

Neurophysiology Department

RHC Glasgow

0141 452 4288 or [paedseeg@ggc.scot.nhs.uk](mailto:paedseeg@ggc.scot.nhs.uk)

## **HOME VIDEO AMBULATORY**

### **GUIDE TO CAMERA SET UP**

- Place camera on a flat surface, for example a table.
- Plug mains cable into wall socket.
- Open the view window (screen) on the side of the camera.
- Press the start/stop button (at back of the camera on the right hand side) REC will be displayed on the screen.
- Please note that the recording will stop if the screen is closed.
- Make sure there is adequate lighting.
- When moving the camera into the bedroom. Make sure it is plugged in and on a flat surface. Press the NIGHTSHOT button on the side of the camera.
- In the morning press the NIGHTSHOT button again.
- The recording will stop if the view window (screen) is closed. If this happens open the view window and press the start/stop button.